

Pediatric Mechanical Ventilation Decision Support User Manual (Lamb Version)

Collaborative Pediatric Critical Care Research Network
Eunice Kennedy Shriver National Institute for Child Health
and Human Development (NICHD)

Julie Beckstrom, M.S.
J. Michael Dean, M.D.
Katherine Sward, Ph.D.

Program Version 0.01
Version Date: August 5, 2012
Printing Date: August 5, 2012

Copyright © 2012. University of Utah School of Medicine on behalf of the Collaborative Pediatric Critical Care Research Network (CPCCRN). All rights reserved.

This user manual describes computer decision support software that is being developed in support of CPCCRN research protocols, and the lead CPCCRN investigators for this research effort are J. Michael Dean, M.D. and Katherine Sward, Ph.D., University of Utah School of Medicine and College of Nursing. The preclinical version of software is being deployed in the Lamb Intensive Care Unit (LICU) in the laboratory of Kurt Albertine, Ph.D., at the University of Utah, in support of NHLBI funded studies on neonatal respiratory failure.

The CPCCRN Clinical Centers are the Children's Hospital of Los Angeles, Children's Hospital of Michigan, Children's Hospital of Philadelphia, Children's Hospital of Pittsburgh, Children's National Medical Center, Phoenix Children's Hospital, and the University of Michigan, and are supported by Cooperative Agreements U10-HD050012, U10-HD050096, U10-HD063108, U10-HD049983, U10-HD049981, U10-HD063114, and U10-HD063106, respectively, from the *Eunice Kennedy Shriver* National Institute for Child Health and Human Development (NICHD).

This document was prepared by the CPCCRN Data Coordinating Center located at the University of Utah School of Medicine, Salt Lake City, Utah. The CPCCRN Data Coordinating Center at the University of Utah is supported by Cooperative Agreement U01-HD049934 from the National Institute for Child Health and Human Development (NICHD). The document was written and typeset using L^AT_EX 2_ε.

THIS PAGE IS INTENTIONALLY BLANK.

Contents

1	Introduction	6
2	Installation	6
3	General Navigation	8
4	Administrative Usage	9
4.1	Creating New Users	9
4.2	Setting Preferences	9
4.2.1	Database Preferences	10
4.2.2	Study Preferences	11
4.3	Database Management	13
4.4	Viewing the Rule Trace	15
5	Normal Usage	16
5.1	Creating, Editing, and Selecting Patients	16
5.2	Entering Laboratory Data	16
5.3	Obtaining Ventilator Decisions	20
6	Summary	23

List of Figures

1	Login screen for decision support program.	7
2	Graphic user interface of decision support software.	8
3	Database menu.	9
4	Dialog window for entering new user information.	10
5	Database preferences dialog.	12
6	Study preferences dialog.	14
7	Study preferences dialog.	15
8	Rule trace view of decisions.	15
9	Entering a new patient into the software.	17
10	List of patients in the database.	18
11	Entering blood gas data.	19
12	Laboratory table view.	19
13	Graphic user interface after pressure control mode selected.	21
14	Graphic user interface for pressure control ventilation.	21

15	Entry screen for ventilator action taken.	22
16	Entry screen when user declines computer recommendation.	24

1 Introduction

The CPCCRN considers numerous scientific concepts for implementation within the network, amongst which have been proposed studies of pediatric mechanical ventilation. The management of mechanical ventilators in the Pediatric Intensive Care Unit (PICU) is complex, both because of the variety of modes of ventilation permitted by current vendors and because of the large number of clinicians involved in the decision making.

[CONSIDER putting in stuff from Alan Morris' work, etc. Kathy might want to write this, or Julie. Idea would be to give about a page of context.]

The software described in this User Manual has been developed in Java, using JBoss Drools as the inference engine that creates advice for the clinician. The software runs on laptop computers and stores the data locally, or optionally may be connected to a central database server. The rules that are used by the inference engine are designed separately from the underlying Java code that provides the user interface.

There are two general modes of use for the software, administrative and bedside user. The purpose of the administrative mode is to allow the creation of user accounts, to enable setup of study preferences, and to manipulate the database settings. This mode also permits visualizing the rules that were fired for specific decisions. The normal bedside user mode prevents accidental destruction of the databases, and does not allow creation of new users.

2 Installation

The software is provided as a zip file that is decompressed with standard software. Your computer must have a Java 1.6 runtime installed or this software will not run. Satisfactory performance has been obtained with vintage 2007 PC laptops. Versions of the software are available for Windows, Linux, or OS X.

Start the application by double-clicking the icon, and a login screen will appear (Figure 1 on the facing page). You will also see (briefly) a window indicating that the knowledge engine is being created, which requires up to a minute. After this disappears, you will be presented with the login dialog, where you enter your name and password. Click on the OK button to start the program.



Figure 1: Login screen for decision support program.

3 General Navigation

When you log in successfully, the software will display an interface similar to Figure 2. This figure shows what an administrative user will see. The normal bedside user view is very similar. The only differences are absence of some of the menus, absence of the Property View, and absence of the Rule Trace tab.




Figure 2: Graphic user interface of decision support software.

4 Administrative Usage

If you login as an administrative user, then you have the ability to create new users, set study preferences, and manipulate the databases. In addition, you are able to view the rule trace, which shows exactly how the decision support was provided.

4.1 Creating New Users

To create a new user, click on the Database menu and select “Add new user” (see Figure 3).





Figure 3: Database menu.

The computer will display a dialog window in which you enter information about the user (Figure 4 on the next page).

You should enter the last and first names of the new user, and then enter the account name which will be used for logging into the software. You will indicate whether the account is an administrative account by checking the Administrative Account checkbox. Enter the password for the user. Note that in the current version of the software, the password cannot be altered and you will need to provide the password to the new user.

4.2 Setting Preferences

If you select the Preferences menu item (location differs by computer platform), you will see that there are four categories of preferences. Two of these are being used by



Add New User to Application

Enter the user information and press the OK button.

Last name:

First name:

Account name:

☐ Administrative Account

Password:

Confirm password:

Cancel OK

Figure 4: Dialog window for entering new user information.

the developers and you should ignore these (Drools Preferences and Pediatric Ventilator Decision Preferences). In the future, the Drools Preferences will control updating decision rules from the Internet, and the Pediatric Ventilator Decision Preferences will allow selection of International Units for international users.

4.2.1 Database Preferences

Select the Database Preferences item from the Preference list. The Database Preferences dialog will be displayed (Figure 5 on page 12). The software has an embedded HSQLDB database for storing data locally on your computer. Optionally, you may choose to store data on a MySQL Server, in which case you should enter the connection settings for your MySQL database. Leave the Hibernate Dialect unchanged. Be careful because the decision support software will not create the MySQL database if it does not exist. In MySQL, you must create an empty database with the name provided in the settings.

The software WILL create the required table structure when you initialize the database (Section 4.3 on page 13).

The target database is a database to which you can transfer data from the main working database. For example, you may store your patient data on the local HSQLDB database or the MySQL database that you set up in the previous step, but wish to send all the data to a central research database (the target). This would be the typical situation in multi-institutional studies implemented by CPCCRN. You would probably set up a MySQL server in your PICU to store all the local data, and the target database would be at the Data Coordinating Center. Encryption options are available to reduce transmission (to the target database) of patient identifiers (Section 4.2.2).

If you are uncertain how to set the database preferences, you may leave them unchanged and the local HSQLDB database will be sufficient. When you are ready to combine data from multiple patients, however, you will need to set up the target database.

After you have finished editing the settings in the dialog, click the Apply button and the OK button. If you wish to restore default settings, you click the Restore Defaults button. Of course, you may Cancel your choices by clicking the Cancel button.

4.2.2 Study Preferences

Select the Study Preferences item from the Preference list. The Study Preferences dialog will be displayed (Figure 6 on page 14). The top section of the dialog determines how study subjects are numbered. The default setting is the simplest and probably the correct setting if you are sharing research data outside your institution. The Prefix of the study site and the Prefix of the study name are used inside the study number, which simplifies combining data from multiple centers.

The software can be run in active or passive mode, determined by the Active Mode checkbox. If this is not checked, then the software is in passive mode. In this mode, users can enter data, but the computer software will not provide any advice. It will ask the user to enter the decision made by the user. This is an important part of the software, because this allows collection of ventilator management data that is being used when there is no computer decision support.

If you check the Active mode checkbox, then the computer program will provide advice when you select the Calculate Decision button, described in Section ?? on page ??.

type filter text

- ▼ Decision Support Preference
 - Database Preferences**
 - Drools Preferences
 - Pediatric Ventilator Decis
 - Study Preferences

Database Preferences

Preference settings for decision support database:

Select database:

☒ Local HSQLDB

☐ MySQL Server

Settings for local HSQLDB database:

JDBC Driver:

Connection URL:

Settings for mySQL database:

JDBC Driver:

Connection URL:

Hibernate Dialect:

Username:

Password:

Settings for TARGET database for transferring data:

JDBC Driver:

Connection URL:

Hibernate Dialect:

Username:

Password:

Restore Defaults Apply

Cancel OK

Figure 5: Database preferences dialog.

Finally, you can turn Encryption On by checking the Encryption On checkbox. When this is checked, then names and medical record numbers on the local database (or local MySQL server) are encrypted with a sophisticated one-way hash prior to transmission to the target database (usually at the Data Coordinating Center). While the same names and medical record numbers will result in the same encryption, it is mathematically impossible to obtain the names or medical record numbers from the encrypted information. The hash is a one way process, so the privacy of your patients is fully protected when this setting is selected.

4.3 Database Management

The Database menu contains several actions in addition to Add new user, which you have already learned how to use. The first one, Initialize Database, will create a new database with the correct structure (tables) for the software to operate. Notice, however, that initializing a database erases its contents (and starts over from scratch). This is permanent (Figure 7 on page 15). The primary reason that this function is not available to normal bedside users is to prevent the accidental destruction of study data with this command.

There is a Add fake patients menu item that is used by the developers to test the software. If you are installing a version of software on a demonstration machine, as you might do to teach people how to use the software, then you can add fake patients with this command. This can allow new users to learn how to use the software without messing up real patient data.

The command to Transfer data to target will move the data on your local database to the target, which is usually at the Data Coordinating Center. This should be self explanatory. It does *not* erase the data from your local machine. If Encryption On has been selected in the Study Preferences, then the patients' names and medical record numbers are encrypted before they are transferred to the target database.

Finally, the administrative user can erase the target database. This is similar to Initialize Database in that it completely erases the previous database. This is a very dangerous command that will probably be removed from the software in later versions. If, for example, the target database at the Data Coordinating Center was connected to your machine, and you selected this option, it would destroy the database at the Data Coordinating Center. This option is present because there are no current studies, and the software is still in development.

The screenshot shows a 'Preferences' window with a sidebar on the left containing a search bar labeled 'type filter text' and a list of preference categories: 'Decision Support Preference', 'Database Preferences', 'Drools Preferences', 'Pediatric Ventilator Decis', and 'Study Preferences' (which is highlighted). The main panel is titled 'Study Preferences' and contains the following sections:

- Preference settings for study information**
 - Select method to create study identification number:**
 - ☒ Create with site and study prefixes
 - ☐ Sequential simple numbers
 - ☐ Free form – user must enter
 - Information below is used to construct record numbers. The prefix strings must be four characters in length.**
 - Last number used:** 0001
 - Name of study site:** University of Utah
 - Prefix of study site (4 characters):** UTAH
 - Name of study:** Demonstration Study
 - Prefix of study name (4 characters):** DEMO
- Status of the decision inference engine.** If active is checked, then the inference engine will fire. If not, the clinician will be asked to enter their actions without receiving advice.
 - ☐ Active Mode
- Set encryption of data exports on or off.** Default is off to facilitate debugging. When on, names and medical record numbers will be encrypted prior to transfer to the MySQL target database.
 - ☐ Encryption On

At the bottom right of the main panel are buttons for 'Restore Defaults' and 'Apply'. At the bottom of the window are buttons for 'Cancel' and 'OK'.

Figure 6: Study preferences dialog.




Figure 7: Study preferences dialog.

4.4 Viewing the Rule Trace

The administrative user can also examine how the software reached its conclusions and constructed its advice to the bedside user, by viewing the Rule Trace (Figure 8). In future releases of this software, we may move this function to the normal bedside user, because there is nothing dangerous about it. The main reason for it to be in the administrative mode is that this may be confusing to users at the bedside.




Figure 8: Rule trace view of decisions.

5 Normal Usage

When you use the software as a normal bedside user, you need to be able to add or select a specific patient, edit patient information, enter laboratory data, and enter information required for each ventilator decision. After the software provides you with advice, then you need to accept or decline the advice. You can also view the previous decisions and all the laboratory results that you have entered into the software.

5.1 Creating, Editing, and Selecting Patients


In the File menu, select New patient. This will bring up a dialog for entering new patient information (Figure 9 on the next page). Notice that the date of birth is automatically initialized with the current date. This will be convenient in the neonatal setting or the Lamb Intensive Care Unit, but in normal PICU usage, you will need to edit and enter the correct birthdate.

When you start the software, there may be no patient selected, or perhaps the incorrect patient is selected. You should click the button to Select (or Change) Active Patient. This button is in the upper left section of the software interface. If there are patients in the database, a listing will be provided (Figure 10 on page 18). Click on the correct patient, and click the OK button. Alternatively, you may double click on the correct patient and it will be selected automatically.

5.2 Entering Laboratory Data

To enter arterial blood gases, you may select the New ABG button on the upper left of the software screen, or you may use the Laboratory menu and select Enter arterial blood gases. Either method results in the Arterial Blood Gas Panel Result dialog (Figure 11 on page 19). The date and time are populated with the current time, and you may need to edit these variables. The required measurements include inspired oxygen as well as the blood gas results.

After laboratory data have been entered, you can view them in the Laboratory Table in the lower section (center) of the screen (Figure 12 on page 19). The serum bicarbonate and base excess or deficit are calculated with standard equations.



The screenshot shows a software window titled "New Research Patients". Below the title bar is a header section with the title "Demographic Information about Study Patient" and a subtitle "Enter demographic information about the patient who is participating in the research study." The main area contains several input fields: "Study ID Number (automatic)" with a value of "DEMOUTAH0001"; "Enter the first name of the patient:" with an empty "First name:" field; "Enter the last name of the patient:" with an empty "Last name:" field; "Enter the date of birth of the patient (MM/DD/YYYY):" with a date of "8/ 4/2012" and a calendar icon; "Enter the medical record number of the patient:" with an empty "Medical record number:" field; "Enter the weight of the patient (kilograms):" with an empty "Weight (kg)::" field; and "Enter the height of the patient (cm):" with an empty "Height (cm)::" field. At the bottom right are "Cancel" and "Finish" buttons.

New Research Patients

Demographic Information about Study Patient

Enter demographic information about the patient who is participating in the research study.

Study ID Number (automatic)

Study ID number: DEMOUTAH0001


Enter the first name of the patient:

First name:

Enter the last name of the patient:

Last name:

Enter the date of birth of the patient (MM/DD/YYYY):

Date of birth: 8/ 4/2012 

Enter the medical record number of the patient:

Medical record number:

Enter the weight of the patient (kilograms):


Weight (kg)::

Enter the height of the patient (cm):

Height (cm)::

Cancel **Finish**


Figure 9: Entering a new patient into the software.



A screenshot of a software window titled "Select the Active Patient". The window contains a table with four columns: "Study ID", "Patient Name", "Medical Record Number", and "Birthdate". The table lists five patients, with the first five rows highlighted in light blue. Below the table are "Cancel" and "OK" buttons.

Study ID	Patient Name	Medical Record Number	Birthdate
ST02CHOM0001	Pediatric, John	12-34-56	December 15, 1999
ST02CHOM0002	Adult, Susan	23-34-56	October 15, 1980
ST02CHOM0003	Pediatric, Steven	34-34-56	November 15, 1999
ST02CHOM0004	Adult, Mary	45-34-56	December 25, 1953
ST02CHOM0005	Adult, Esther	56-34-56	December 15, 1969

Figure 10: List of patients in the database.



Arterial Blood Gas Panel Result

Enter the date, time, and values for the arterial blood gas.

Date: 8/ 5/2012

Time: 5:16 PM

FiO2: (21 to 100)


pH:

pCO2: mmHg

pO2: mmHg

Cancel OK

Figure 11: Entering blood gas data.



Decision Table		Laboratory Table							
Laboratory date		FiO2	pH	pCO2	pO2	PaO2/FiO2	HCO3	BE/BD	
<input type="checkbox"/>	Sat Aug 04 09:27:51 MDT 2012	35	7.23	44	32	91	17.8	-8.5	

Figure 12: Laboratory table view.

5.3 Obtaining Ventilator Decisions

The purpose of the software is to handle ventilator decisions. The first step is to select the mode of ventilation. There are five modes supported by this software:

- Pressure Control
- Pressure Regulated Volume Control (PRVC)
- High Frequency Oscillatory Ventilation (HFOV)
- Volume Control
- Extubated

The mode can be selected by clicking on one of the radio buttons, or by selecting the appropriate tab. In the Lamb ICU, you should select the Pressure Control mode. Now the whole interface will be filled in, similar to [Figure 13 on the next page](#).

Let's look at the variables for Pressure Control in more detail ([Figure 14 on the facing page](#)).

At the top, there is a button for Refresh Date Time, as well as separate controls that you can use to edit the date and time. The Refresh Date Time button will set these controls to the current date and time. On the upper right, the New ABG will allow entry of a blood gas.

In the middle section, there is set of tab folders, one for each mode of ventilation. The most recent blood gas is displayed; if there is no previous blood gas then these fields will be empty. The subsequent fields that are in yellow are required. After you have entered the required fields, then the Calculate Decision button will become enabled. When that happens, then you can click on that button to obtain advice for your patient.

NOTE: The current software will be changed so that SpO2 is optional (since it is not available in the Lamb ICU). In addition, we know that we need to add an inspiratory time field. This will be done before the software is officially deployed in the LICU.

If the software is in passive mode, then when you click the Calculate Decision button, you will be asked to enter the actual ventilator decision that you chose to make ([Figure 15 on page 22](#)). Filling this in accurately is important because the software development team is using your current ventilator management to help define rules to tell the computer how to develop advice for future users.

[illegible]

Figure 13: Graphic user interface after pressure control mode selected.

Refresh Date Time

Observation Date: 8/ 4/2012

Time (hh:mm): 9:25:01 AM

New ABG

☒ Pressure Control

☐ PRVC

☐ HFOV

☐ Volume Control

☐ Extubated

Pressure Control
 Pressure Regulated VC
 High Frequency Oscillatory
 Volume Control
 Extubated

pH: Not available
 PaO2: Not available
 pCO2: Not available

Not available
 Not available
 Not available

FiO2:
 SpO2:
 ETCO2:

Ventilator Rate:
 VT (total ml):
 Calculated ml/kg

PEEP:
 PIP:
 MAP:

Figure 14: Graphic user interface for pressure control ventilation.

The screenshot shows a software window titled "Clinical actions taken at this time." with a subtitle "Please enter the decisions that were made by the clinical team in the fields below. You may also add descriptive comments." Below the subtitle are two radio buttons: "No changes at this time" and "Changes as shown below:". Under "Changes as shown below:" is a section titled "Ventilator Changes" containing eight checkboxes, each followed by a parameter name and a text input field. The first checkbox is for "FiO2 (%)" with a "Current value: 21%" displayed to its right. The other parameters are Ventilator rate, PEEP, PIP, Tidal volume (total cc), Amplitude (HFOV only), Frequency (HFOV only), and Mean airway pressure. Below the "Ventilator Changes" section is a section titled "Descriptive comments (optional)" with a large text area. At the bottom of the window, there is a footer with version information and an "OK" button.

Clinical actions taken at this time.

Please enter the decisions that were made by the clinical team in the fields below. You may also add descriptive comments.

☐ No changes at this time

☐ Changes as shown below:

Ventilator Changes

☐ FiO2 (%): Current value: 21%

☐ Ventilator rate:

☐ PEEP:

☐ PIP:

☐ Tidal volume (total cc):

☐ Amplitude (HFOV only):

☐ Frequency (HFOV only):

☐ Mean airway pressure:

Descriptive comments (optional)

Ventilator Decision Support User Manual Version 0.01.01
Manual Version Date: August 5, 2012
The NICHD Collaborative Pediatric Critical Care Research Network

OK

Figure 15: Entry screen for ventilator action taken.

If the software is in active mode, then when you click the Calculate Decision button, you will be presented with the Decision Recommended by Protocol in the middle panel, and there should be an explanation in the panel below it. At this point, you can Accept Decision by clicking the appropriate button, in which case the decision is stored. At this point, the clock in the lower left of the screen will be reset and will count down to the time that you should enter the next ventilator decision (or status).

If you do not agree with all of the advice, you should click the Decline Decision button, and then you will be presented with a dialog to ask you what you did, and why you disagree with the computer (Figure [16 on the following page](#)). It is important for you to realize that you are almost certainly correct if you do not agree with the computer. The information you enter at this point is critical for the developers to improve the software so that it makes better decisions.

6 Summary

This computer decision support software is a computer program in evolution. You should not hesitate to contact the authors with suggestions or questions. Our immediate goals are to deploy working software in the Lamb ICU at the University of Utah so that we can perfect the usability of the software, as well as develop rules that will enable the software to provide reasonable recommendations for ventilator management. In the future, this software will be used in human subjects in the Pediatric ICU setting. Your help now will enable this software to reduce morbidity and mortality from ventilatory failure in children.

The screenshot shows a software window titled "Ventilator Decision Support User Manual Version 0.01.01". The window has a title bar with three standard OS window controls (minimize, maximize, close). The main content area is divided into several sections. At the top, a header section with a blue gradient background contains the text "Reason for choosing to decline the recommended decision." followed by "Please describe your reasoning so that we can improve the software." Below this, there are two radio buttons: "No changes at this time" and "Changes as shown below:". The "Changes as shown below:" option is selected. Underneath, a section titled "Ventilator Changes" contains a list of eight parameters, each with a checkbox and a text input field. The parameters are: FiO2 (%), Ventilator rate, PEEP, PIP, Tidal volume (total cc), Amplitude (HFOV only), Frequency (HFOV only), and Mean airway pressure. To the right of the FiO2 input field, the text "Current value: 21%" is displayed. Below the "Ventilator Changes" section, there is a section titled "Reasons for declining decision (REQUIRED):" followed by a large, empty rectangular text area for user input. At the bottom of the window, there is a footer section containing the text "Ventilator Decision Support User Manual Version 0.01.01", "Manual Version Date: August 5, 2012", and "The NICHD Collaborative Pediatric Critical Care Research Network". To the right of the footer text is an "OK" button.

Reason for choosing to decline the recommended decision.
Please describe your reasoning so that we can improve the software.

☐ No changes at this time
☒ Changes as shown below:

Ventilator Changes

☐ FiO2 (%): Current value: 21%
☐ Ventilator rate:
☐ PEEP:
☐ PIP:
☐ Tidal volume (total cc):
☐ Amplitude (HFOV only):
☐ Frequency (HFOV only):
☐ Mean airway pressure:

Reasons for declining decision (REQUIRED):

Ventilator Decision Support User Manual Version 0.01.01
Manual Version Date: August 5, 2012
The NICHD Collaborative Pediatric Critical Care Research Network

OK

Figure 16: Entry screen when user declines computer recommendation.